

pH and ORP Calibration Buffer Sachets



No pre-mixed Buffer degradation

Mix with 100 mL of DI water for instant fresh standard
pH 4, pH 7 and pH 10 sachets are standard

ORP 90mV and 255 mV sachets are standard

Save on transportation costs compared to liquid

These pH and ORP buffer sachet pillows offer an excellent alternative to pre-mixed buffer calibration standards. Avoid the expensive transport of liquid buffer. These buffer sachets provide a ready to mix 100 mL sample of buffer calibration solutions. Simply mix the contents of the buffer pillow with 100 mL of distilled water to prepare an instant calibration sample.

The standard sachet pH buffers offered are 4, 7 and 10. For ORP (redox), a 90mV and 255mV sachet pillows are available.

Special buffer formulations are available based on required value and volume.

Ordering Information

Buffer Calibration Sachet pillows: (Mix with 100mL DI water)

Part Number	Description
BP104-100	pH Buffer 4.01 Buffer sachet, dry powder sealed in foil package,. Rehydrate to 100mL.
BP107-100	pH Buffer 7.00 Buffer sachet, dry powder sealed in foil package,. Rehydrate to 100mL.
BP110-100	pH Buffer 10.00 Buffer sachet, dry powder sealed in foil package,. Rehydrate to 100mL.
BP107QH-100	ORP 90mV Standard sachet, dry powder. Rehydrate to 100mL.
BP104QH-100	ORP 255mV Standard sachet, dry powder. Rehydrate to 100mL.

NIST Traceable